

A guided tour for analyzing and quantifying single- molecule localization microscopy data

Part 1: Analysing raw SMLM data

NEUBIAS@Home webinar


11/06/2020

Siân Culley a.k.a. @SuperResoluSian

Overview



- Basic principles of SMLM analysis
- Analysis for sparse datasets
- Analysis for dense datasets
- Quality control
- 3D localisation

Basic principles of SMLM



Detection and localization

- Approximate molecule locations = detection
- Accurate determination of molecule centre = localization



(x, y) coordinates of molecule within image

e.g.
(45.29nm, 53.75nm)

ThunderSTORM

- Fiji plugin, all GUI-based
- Can handle 2D and 3D data
- Has options for high density analysis
- Easy to export particles table
- Generates protocol .txt files

nature methods

ANALYSIS

<https://doi.org/10.1038/s41592-019-0364-4>

Corrected: Publisher Correction

Super-resolution fight club: assessment of 2D and 3D single-molecule localization microscopy software

Daniel Sage^{1,22*}, Thanh-An Pham^{1,22}, Hazen Babcock², Tomas Lukes^{3,4}, Thomas Pengo⁵, Jerry Chao^{6,7}, Ramraj Velmurugan^{7,8}, Alex Herbert⁹, Anurag Agrawal¹⁰, Silvia Colabrese¹¹, Ann Wheeler¹², Anna Archetti¹³, Bernd Rieger¹⁴, Raimund Ober^{6,7,15}, Guy M. Hagen¹⁶, Jean-Baptiste Sibarita^{17,18}, Jonas Ries¹⁹, Ricardo Henriques²⁰, Michael Unser¹ and Seamus Holden^{21,22*}

Sage et al, Nature Methods (2019)

ThunderSTORM: a comprehensive ImageJ plug-in for PALM and STORM data analysis and super-resolution imaging

Martin Ovesný, Pavel Křížek, Josef Borkovec, Zdeněk Švindrych and Guy M. Hagen*

Institute of Cellular Biology and Pathology, First Faculty of Medicine, Charles University in Prague, Prague 12800, Czech Republic


Associate Editor: Jonathan Wren

Supplementary data – bonus maths, user manual

<https://github.com/zitmen/thunderstorm/wiki>

Fiji install:

Help > Update... > Manage Update Sites




Hohlbein Lab
update site

Data 'density'

Real experimental datasets from <http://bigwww.epfl.ch/smlm/datasets/index.html>

Low density
(sparse)


- Emitters are all very well separated
- Few emitters per μm^2
- Typically acquired over long period of time (i.e. large number of frames)
- Drift can be severe



2 μm

High density

- Emitters overlap
- Several emitters per μm^2
- Often fewer frames
- From live cells? Watch out for movement



2 μm

'emitter' = fluorophore in 'on' state

Particles tables

ThunderSTORM: results

id	frame	x [nm]	y [nm]	sigma ...	intensity [pho...	offset [photon]	bkgstd [photon]	chi2	uncertainty_...
1	1	417.468	2640.243	223.509	10258.275	44.722	46.445	3323.058	16.741
2	1	481.68	7056.711	159.844	5231.573	80.517	79.321	9647.297	28.435
3	1	506.473	1580.226	153.191	7259.846	80.155	40.213	3105.818	10.029
4	1	669.437	6352.407	128.805	1669.543	78.791	19.009	1276.271	15.073
5	1	624.57	7281.365	168.006	5891.138	59.687	70.991	7192.431	24.982
6	1	792.393	547.101	152.966	4676.928	125.514	67.415	7679.431	24.901
7	1	926.555	5699.312	162.931	3687.128	71.528	28.936	2198.012	15.912
8	1	1033.421	2671.934	289.785	8459.08	35.94	24.049	1469.915	18.011
9	1	1094.673	2843.861	145.309	1351.337	119.062	44.481	3936.051	51.264
10	1	1372.559	3299.007	170.209	7062.607	70.419	15.525	523.907	5.881
11	1	1348.05	11233.328	155.389	2743.277	58.665	27.535	1996.595	18.464
12	1	1600.821	1641.058	150.626	2556.421	67.228	17.356	820.734	12.487
13	1	1656.894	1201.632	172.935	5218.995	68.678	13.517	483.853	7.076
14	1	1622.658	8669.813	133.462	2179.423	77.529	35.173	2802.113	21.755
15	1	1742.882	2685.776	139.385	2338.193	82.256	10.663	337.574	8.315
16	1	1723.45	5215.905	143.719	1401.231	63.303	7.832	277.202	10.909
17	1	1945.564	6888.527	178.277	3650.194	61.406	37.683	3457.906	24.342
18	1	1846.094	12430.179	147.388	2036.822	61.811	22.148	1310.547	18.259
19	1	1923.11	9138.857	133.476	5498.35	71.501	19.218	924.297	5.621
20	1	2189.383	7182.671	172.179	4087.851	53.205	37.227	3116.562	20.171
21	1	2343.263	3617.819	119.901	709.139	68.663	5.785	154.815	11.86
22	1	2443.59	313.595	143.091	4378.591	95.276	38.468	3472.473	13.741
23	1	2531.572	2171.223	198.974	3405.234	60.527	29.085	2358.366	25.157
24	1	2633.573	2235.262	142.937	1772.464	86.638	36.36	3707.549	31.268
25	1	2627.909	4278.172	136.752	3577.781	70.637	17.783	982.555	7.992
26	1	2813.703	10091.223	141.82	4869.598	58.596	20.503	1092.351	7.221
27	1	2745.725	12256.168	122.87	1646.154	68.07	10.394	475.463	8.928
28	1	2788.333	4687.385	175.656	2821.404	49.672	31.47	2321.055	25.637
29	1	2894.192	11634.63	137.295	5426.208	63.545	24.844	1432.882	7.189
30	1	2920.342	2761.895	126.548	4446.844	74.274	15.003	727.497	5.248
31	1	2962.616	691.525	153.818	2965.512	55.643	19.615	976.249	12.52
32	1	3009.069	8109.174	117.923	1277.683	82.249	21.674	1522.845	18.522
33	1	3028.664	8574.056	123.54	2865.356	80.048	17.519	978.756	8.083
34	1	3175.319	7603.575	157.209	2051.041	65.701	22.3	1772.093	20.562

Filter: Density filter Remove duplicates Merging Drift correction Z-stage offset

Filter: Apply Restrict to ROI



Post-processing history: -

☒ Preview

central coordinate

width of Gaussian

amplitude
and
offset of
Gaussian



χ^2 goodness-of-fit

uncertainty

$$\langle (\Delta x)^2 \rangle = \frac{2\sigma^2 + a^2/12}{N} + \frac{8\pi\sigma^4 b^2}{a^2 N^2}$$

Analysis woes for high-density data...





Raw data

Detections





Rendered image

Widths of fitted Gaussians
(‘Plot Histogram’ button)

Low
density






High
density



High density data in ThunderSTORM

- Just filtering out 'bad' detections?
- ThunderSTORM: multi-emitter fitting



Batch processing with ThunderSTORM

<https://jalink-lab.github.io/>

Jalink-lab.github.io

SMLM Fiji plugins and macros

Single-Molecule Localization Microscopy Tools

Here you find several Fiji scripts and plugins to perform automated (batch) processing of Single-Molecule Localization Microscopy (SMLM) datasets, using the ImageJ plugin ThunderSTORM. ThunderSTORM is a very useful tool for analysis and visualization of localization microscopy data. However, it doesn't have much functionality for batch processing, and it misses other necessary processing steps, like temporal background subtraction and chromatic aberration correction. We have developed some tools to automate these processes.


Download the complete package:

1. ImageJ1 macro SMLM_process_folder.ijm
2. Plugin Temporal Median Background Subtraction
3. Plugin Chromatic Aberration Correction
4. Plugin ImageJSON

Installation instructions

ThunderSTORM alternatives...


- HD algorithms...




Pre-processing high density data

- HAWK wavelet filtering
- Makes datasets more 'sparse'

Example raw high-density frame



Example high-density frame after HAWK





nature | **methods**

BRIEF COMMUNICATION

<https://doi.org/10.1038/s41592-018-0072-5>

Artifact-free high-density localization microscopy analysis

Richard J. Marsh, Karin Pfisterer, Pauline Bennett, Liisa M. Hirvonen , Mathias Gautel,
Gareth E. Jones and Susan Cox *

Download from: www.coxphysics.com

Non-particles table approaches

- SOFI, SRRF...

RESEARCH ARTICLE

Fast, background-free, 3D super-resolution optical fluctuation imaging (SOFI)

T. Dertinger, R. Colyer, G. Iyer, S. Weiss, and J. Enderlein

PNAS December 29, 2009 106 (52) 22287-22292; <https://doi.org/10.1073/pnas.0907866106>

Edited by John W. Sedat, University of California, San Francisco, CA, and approved October 29, 2009 (received for review July 15, 2009)

Open Access | Published: 12 August 2016

Fast live-cell conventional fluorophore nanoscopy with ImageJ through super-resolution radial fluctuations


Nils Gustafsson, Siân Culley, George Ashdown, Dylan M. Owen, Pedro Matos Pereira & Ricardo Henriques ✉

Nature Communications **7**, Article number: 12471 (2016) | [Cite this article](#)

10k Accesses | 141 Citations | 95 Altmetric | [Metrics](#)



SRRF analysis of high density data




- SRRF is available as a Fiji plugin
- Used for very high density data
- Image \approx probability map, not intensity or molecule coordinates


Artefacts in SMLM reconstructions

- Merging of closely-separated structures
- Incomplete representation of structures
- Intensity non-linearities


Simulated raw data



SMLM reconstruction




Real structure



2μm


Low density data –
filtering around
fiducial marker



2μm

Using SQUIRREL to assess image quality

- Can compare different algorithm results
- Can map resolution across image
- Can test how many frames you need to represent an image



Published: 19 February 2018

Quantitative mapping and minimization of super-resolution optical imaging artifacts

Siân Culley, David Albrecht, Caron Jacobs, Pedro Matos Pereira, Christophe Leterrier

✉, Jason Mercer ✉ & Ricardo Henriques ✉

Nature Methods **15**, 263–266(2018) | [Cite this article](#)

3013 Accesses | **64** Citations | **160** Altmetric | [Metrics](#)

NanoJ-Core,
NanoJ-SQUIRREL
update sites

Manage update sites		
...	Name	URL
<input type="checkbox"/>	Molography	https://sites.imagej.net/Volker
<input type="checkbox"/>	MoMA	https://sites.imagej.net/MoMA
<input type="checkbox"/>	Morphology	https://sites.imagej.net/Landir
<input checked="" type="checkbox"/>	MOSAIC ToolSuite	https://mosaic.mpi-cbg.de/Doi
<input type="checkbox"/>	MPIBPC	https://sites.imagej.net/MPIBPC
<input type="checkbox"/>	MS-ECS-2D	https://sites.imagej.net/MS-ECS-2D
<input type="checkbox"/>	MTrack	https://sites.imagej.net/MTrack
<input type="checkbox"/>	Multifrac	https://sites.imagej.net/Multifrac
<input type="checkbox"/>	Multi-Template-Matching	https://sites.imagej.net/Multi-Template-Matching
<input type="checkbox"/>	N5	https://sites.imagej.net/N5/
<input checked="" type="checkbox"/>	NanoJ-Core	https://sites.imagej.net/NanoJ-Core
<input checked="" type="checkbox"/>	NanoJ-SQUIRREL	https://sites.imagej.net/NanoJ-SQUIRREL
<input type="checkbox"/>	NanoJ-SRRF	https://sites.imagej.net/NanoJ-SRRF
<input type="checkbox"/>	NanoJ-VirusMapper	https://sites.imagej.net/NanoJ-VirusMapper
<input type="checkbox"/>	NEUBIAS	https://sites.imagej.net/NEUBIAS
<input type="checkbox"/>	Neuroanatomy	https://sites.imagej.net/Neuroanatomy
<input type="checkbox"/>	NeuroCyto Lab	https://sites.imagej.net/NeuroCyto Lab

Latest release

v1.1-alpha
12e40f9

Compare

GPU-independent error mapping

superresolusian released this on Jul 12, 2019

Test release for users whose graphics cards object to NanoJ libraries.

Assets 3

NanoJ_SQUIRREL.jar


Source code (zip)

Source code (tar.gz)


<https://github.com/superresolusian/mypluginsbaby/releases/>

3D data

- You can't get axial localizations unless you used specific optics to encode axial information
- You need to have calibration data for your microscope



Algorithms with 3D data



- ThunderSTORM is good for astigmatism
- SMAP can do all 3D modalities (very high performance for biplane)


Published: 09 April 2018

Real-time 3D single-molecule localization using experimental point spread functions

Yiming Li, Markus Mund, Philipp Hoess, Joran Deschamps, Ulf Matti, Bianca Nijmeijer, Vilma Jimenez Sabinina, Jan Ellenberg, Ingmar Schoen & Jonas Ries

Nature Methods **15**, 367–369(2018) | [Cite this article](#)

3400 Accesses | 36 Citations | 55 Altmetric | [Metrics](#)



General advice for SMLM analysis

- Walk before you run
- Test analysis on small crops of data before running on whole datasets
- Get familiar with one piece of software

Report what you did!

- What algorithm did you use?
- What parameters did you use?
- How are your images rendered?
- Did you test any other algorithms/parameters – how did you decide which was best?